

# P767 Quantum Blue® Infliximab POC User Performance Evaluation

Benoit Fellay<sup>a</sup>, Maria Kopertowska-Majchrzak<sup>b</sup>, Emanuel Burri<sup>c</sup>, Marie-Eve Ueberschlag<sup>d</sup>, Alicja Ritz<sup>d</sup>, Thomas Schuster<sup>d</sup>

<sup>a</sup> Hôpital Cantonal Fribourg, Fribourg, Switzerland

<sup>b</sup> Wielospecjalistyczny Szpital Wojewódzki, Gorzów Wlkp, Poland

<sup>c</sup> Kantonsspital Baselland, Liestal, Switzerland;

<sup>d</sup> BÜHLMANN Laboratories AG, Schoenenbuch, Switzerland ([ts@buhlmannlabs.ch](mailto:ts@buhlmannlabs.ch))



## BACKGROUND

The objective of the user performance evaluation was to demonstrate the ease-of-use of the Quantum Blue® Infliximab POC system, consisting of the Quantum Blue® Infliximab test and the Dilution Set, to allow non-laboratory professionals to independently and correctly determine infliximab concentrations starting from an existing patients' serum sample within 15 minutes under actual conditions of use and without further laboratory equipment.

## METHODS

The ability of point-of-care users to obtain correct results was evaluated by testing result agreement between users' Quantum Blue® Infliximab POC measurements, obtained using three kit lots, in total, and laboratory reference values (RIDASCREEN® IFX Monitoring ELISA (R-Biopharm)) for a set of 40 clinical serum samples (Fig. 1). To demonstrate the ease-of-use of the test for non-laboratory professionals, the performance of the users was compared to that of laboratory personnel at BÜHLMANN who performed the same measurements with the Quantum Blue® Infliximab POC system. Test robustness, ease and comfort of use as well as the clearness of the given instructions was further assessed in a questionnaire. Three POC sites, in three geographically distinct locations participated in this study. Operators were non-laboratory medical personnel such as nurses, medical practice assistants or physicians. Two operators were recruited per site. Sites: 1) Hôpital Cantonal Fribourg, Fribourg, Switzerland, 2) Wielospecjalistyczny Szpital Wojewódzki, Gorzów Wlkp, Poland, and 3) Kantonsspital Baselland, Liestal, Switzerland.

## RESULTS

None of the six non-laboratory professionals of the three POC sites, received a false-positive or false-negative result, based on an optimal therapeutic window of 3 to 7 µg/mL (Vande Casteele et al., 2015). Overall, non-laboratory professionals at the POC sites received comparable results as the laboratory professionals at BÜHLMANN (Fig. 2 and 3). Bias at 3 and 7 µg/mL, clinical decision points for therapeutic drug monitoring, when compared to laboratory reference values, were determined to 4.8 % and 7.4 % (site 1, Fig. 3), 2.4 % and 5.8 % (site 2) as well as to 12.9 % and 17.0 % (site 3). The total agreement of non-laboratory professionals' results with reference infliximab values was 82.3 % (site 1, Fig. 5), 80.8 % (site 2, Fig. 6) and 83.8 % (site 3, Fig. 7) and comparable between sites (Fig. 4). Overall the non-laboratory professionals' assessment of the POC assay in terms of the robustness, ease and comfort of use was very positive.




Fig. 1: Study set-up




Fig. 2: Passing-Bablok regression analysis of professional users (setting and lots as site 1, slope of 0.97)




Fig. 3: Passing-Bablok regression analysis for non-laboratory personnel (site 1, slope of 1.09)




Fig. 4: Bland-Altman analysis (site 1) revealed a bias over the measuring range of 8.58%.

Test	Category	Comparative Method		
		low	optimal	high
Quantum Blue® Infliximab	low	30.4%	1.3%	0.0%
	optimal	7.6%	20.3%	2.5%
	high	0.0%	6.3%	31.6%

82.3%

Fig. 5: Diagnostic agreement for non-laboratory personnel (site 1 vs. comparative method)

Test	Category	Comparative Method		
		low	optimal	high
Quantum Blue® Infliximab	low	26.9%	3.8%	0.0%
	optimal	10.3%	20.5%	1.3%
	high	0.0%	3.8%	33.3%

80.8%

Fig. 6: Diagnostic agreement for non-laboratory personnel (site 2 vs. comparative method)

Test	Category	Comparative Method		
		low	optimal	high
Quantum Blue® Infliximab	low	28.8%	2.5%	0.0%
	optimal	8.8%	22.5%	2.5%
	high	0.0%	2.5%	32.5%

83.8%

Fig. 7: Diagnostic agreement for non-laboratory personnel (site 3 vs. comparative method)

## CONCLUSIONS

The outcome of this study suggests that the Quantum Blue® Infliximab POC System, which determines infliximab levels in serum specimens, is easy-to-use, the given instructions are comprehensive, and the results are comparable between different POC sites as well as between POC sites and laboratories.