

BioTite™ Containers

Leak-proof Containers

Stackable BioTite™ containers are ideal for collection and transport of fluids or solids. The dual-threaded, polyethylene caps provide a super-positive seal with the polypropylene jars, tested at up to 500kPa, BioTites are ideal for use in pneumatic transport systems.

Tough and shatter-resistant with easy to read volume graduation the BioTites design also incorporates a reinforcement ring to strengthen the vial side-wall guarding against squeeze related leakage.

BioTites are available sterile, with tamper evident label or non-sterile, without label.

- Stackable, therefore suitable for automation
- Dual-threaded caps for a leak-resistant seal
- Easy-grip top and bottom for safe handling, even with gloves
- Super-positive seal engages in last 12 degrees of turn
- Reinforcement ring strengthens the vial side-wall, thus preventing squeeze related leakage
- Easy-to-read volume graduations

BioTite™ Leak-proof Containers

Catalogue Number	Description	Pack Size
------------------	-------------	-----------

Sterile with Tamper-evident Label, Pre-capped

LW5465	60ml	500
LW5475	90ml	400
LW5495	120ml	300

Non-Sterile, No Label, Separate Caps

LW5450	20ml	1000
LW5455	40ml	600
LW5460	60ml	500
LW5470	90ml	400
LW5490	120ml	300

