

Proline® Plus – the next generation

- Classic Styling
- Innovative Technology

The Biohit Proline pipette has long been a favourite in laboratories worldwide. Thousands of users have appreciated its quality and durability, giving them years of reliable service.

You'll find the same high standards built into the new Proline Plus pipette. Behind its smart, classic design lies state of the art technology with an emphasis on robustness, ergonomics and efficiency in liquid handling. Available in a full range of single and multichannel versions, Proline Plus is ideal for use in routine labs as well as research applications.



Volume setting with click stops

Suitable for both left and right-handed users

Thermal insulation of the internal components ensures hand warmth does not affect pipetting results

Colour coding on display corresponds to the colour coding of Biohit tip trays

Good chemical and UV resistance

Fits most commonly used lab tubes, thanks to the slim ejector collar

Safe-cone Filters to prevent contamination

Robust design



Fully Autoclavable



CE

IVD

Unique Optoload function on multichannels for easy tip loading and ejection

Special Price Offer

Experience the benefits for yourself.

Take advantage of our fantastic special offer prices on all Proline Plus variable volume pipettes.

Up to 53% discount off list prices!

See overleaf for more details

Proline® Plus Single and Multichannel pipettes

Performance Specifications and Ordering Table

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Vol (µl)	Inacc. (%)	Impr. (%)	List Price (£)	Special Offer Price (£)
728010	1	0.1-3 µl	0.002	3 1.5 0.3	1.30% 2.40% 10.00%	0.80% 1.60% 6.00%	170.50	80.00
728020	1	0.5-10 µl	0.01	10 5 1 0.5	1.00% 1.50% 2.50% 5.00%	0.60% 1.00% 1.50% 4.00%	170.50	80.00
728030	1	2-20 µl	0.02	20 10 2	0.90% 1.20% 3.00%	0.40% 1.00% 2.00%	162.25	80.00
728050	1	10-100 µl	0.1	100 50 10	0.80% 1.00% 2.00%	0.15% 0.40% 1.00%	162.25	80.00
728060	1	20-200 µl	0.2	200 100 20	0.60% 0.80% 2.00%	0.15% 0.30% 0.80%	162.25	80.00
728070	1	100-1000 µl	1	1000 500 100	0.60% 0.60% 1.00%	0.20% 0.20% 0.40%	162.25	80.00
728080	1	500-5000 µl	10	5000 2500 500	0.50% 0.60% 2.00%	0.20% 0.30% 0.60%	162.25	80.00
728090	1	1-10 ml	20	10000 5000 1000	0.60% 1.20% 3.00%	0.20% 0.30% 0.60%	170.50	80.00
728120	8	0.5-10 µl	0.01	10 5 1	1.50% 2.50% 4.00%	1.00% 2.50% 4.00%	459.80	233.00
728130	8	10-100 µl	0.1	100 50 10	0.70% 1.00% 3.00%	0.25% 0.70% 1.50%	459.80	233.00
728140	8	30-300 µl	0.2	300 150 30	0.60% 1.00% 2.00%	0.25% 0.50% 1.00%	459.80	233.00
728220	12	0.5-10 µl	0.01	10 5 1	1.50% 2.50% 4.00%	1.00% 2.50% 4.00%	554.40	295.00
728230	12	10-100 µl	0.1	100 50 10	0.70% 1.00% 3.00%	0.25% 0.70% 1.50%	554.40	295.00
728240	12	30-300 µl	0.2	300 150 30	0.60% 1.00% 2.00%	0.25% 0.50% 1.00%	554.40	295.00

Please quote promotional code QM887A when placing your order.

For more information or to speak to an Alpha representative, call us free on 0800 387732.

Time limited offer, not valid in conjunction with any other offer or quotation. Terms and conditions apply.