

Alpha Pipette Servicing and Calibration

Alpha Laboratories is highly regarded as a supplier of liquid handling products including quality pipettes, tips and other consumables. To compliment the range, a full set of service and calibration options is available covering all popular brands of pipettes and dispensers.

Regular service and calibration will ensure that the precision and accuracy of pipettes are maintained and your laboratory's results are not compromised. In partnership with an experienced calibration team, Alpha offers different levels of service from a simple performance check to a full UKAS calibration with appropriate certification. They can be performed either on-site or at our workshop depending on your requirements.



Workshop Service

Send your pipettes into our workshop where qualified technicians will perform the service and calibration in a purpose-built, environmentally controlled area. Prices for the various options are shown in the price guide overleaf. Simply choose the one you require, download and complete a workshop service request form at www.alphalabs.co.uk/pipette-service and send it to Alpha Pipette Service, (address details are on the form) together with your decontaminated pipettes. We recommend that you use a traceable delivery service to ensure that your pipette reaches us safely. If your pipettes require parts or repairs which are extra to the price of your chosen option, you will be contacted for authorisation to proceed.



On-Site Service

Our technicians will bring their equipment to you and perform an on-site pipette clinic. This is ideally suited to larger quantities of pipettes and can be set up as a one-off visit or on a scheduled basis. The price guide overleaf gives an indication of the charging structure but please contact us to discuss your needs and obtain a customised quotation.

Calibration Options

QuikCheck

A quick and simple performance check for when extensive measurements are not required. An ISO calibration certificate is included.

FullCheck

This service includes cleaning, servicing and a full range performance check. An ISO calibration certificate is supplied. The costs of any spare parts are charged separately.

UKASCal

Accreditation by UKAS gives you the confidence of full pipette traceability. Service level is the same as FullCheck but issued with a UKAS calibration certificate.

FixCheck

This service is for pipettes that are no longer performing within the manufacturer's tolerances or are damaged. It includes dismantling, cleaning, servicing and a full range performance check. An ISO calibration certificate is supplied. All spare parts are detailed and charged separately. This service is not available at on-site clinics.

Further Information

If you require more details on the calibration options or a quotation for on-site pipette clinics, please contact us on 0800 387732 or e-mail: pipette-service@alphalabs.co.uk

PRICE GUIDE

We have a convenient and cost-effective charging structure which is based on the quantity and type of pipettes, the calibration option required and where the work is to be carried out. Workshop prices include return carriage to a UK address. On-site prices include travel and labour and apply to a minimum of 25 pipettes.



	No. of Readings	Test Points	1-50		51-100		101+	
			Workshop	On-site*	Workshop	On-site	Workshop	On-site
Fixed Volume Single Channel	4	4 @ Nominal	£16.20	£16.00	£14.25	£14.00	£13.75	£12.00
Variable Volume Single Channel	4	4 @ Mid Point	£16.20	£16.00	£14.25	£14.00	£13.75	£12.00
Dispenser	4	4 @ Mid Point	£16.20	£16.00	£14.25	£14.00	£13.75	£12.00
Multi 8 Channel	32	4/chan @ Mid Point	£29.75	£28.00	£26.25	£25.00	£24.75	£21.00
Multi 12 Channel	48	4/chan @ Mid Point	£43.75	£33.00	£37.25	£29.00	£36.25	£27.00



	No. of Readings	Test Points	1-50		51-100		101+	
			Workshop	On-site*	Workshop	On-site	Workshop	On-site
Fixed Volume Single Channel	10	10 @ Nominal	£22.25	£17.00	£19.75	£16.00	£19.25	£15.00
Variable Volume Single Channel	18	6 @ Min/Mid/Max	£31.25	£26.00	£27.25	£24.00	£26.25	£22.50
Dispenser	18	6 @ Min/Mid/Max	£31.25	£26.00	£27.25	£24.00	£26.25	£22.50
Multi 8 Channel	144	6/chan @ Min/Mid/Max	£44.25	£33.00	£38.25	£29.00	£37.25	£27.50
Multi 12 Channel	216	6/chan @ Min/Mid/Max	£66.25	£36.00	£57.25	£32.00	£55.25	£30.00



	No. of Readings	Test Points	1-50		51-100		101+	
			Workshop	On-site*	Workshop	On-site	Workshop	On-site
Fixed Volume Single Channel	10	10 @ Nominal	£32.25	£26.00	£27.75	£24.00	£26.75	£22.00
Variable Volume Single Channel	30	10 @ Min/Mid/Max	£52.25	£30.00	£44.75	£28.00	£43.75	£26.00
Dispenser	30	10 @ Min/Mid/Max	£52.25	£30.00	£44.75	£28.00	£43.75	£26.00
Multi 8 Channel	240	10/chan @ Min/Mid/Max	£80.00	£80.00	£69.25	£76.00	£67.25	£74.00
Multi 12 Channel	360	10/chan @ Min/Mid/Max	£121.25	£100.00	£103.25	£95.00	£98.75	£90.00



	No. of Readings	Test Points	Workshop
Fixed Volume Single Channel	10	10 @ Nominal	£33.25
Variable Volume Single Channel	30	10 @ Min/Mid/Max	£55.25
Dispenser	30	10 @ Min/Mid/Max	£55.25
Multi 8 Channel	240	10 /chan @ Min/Mid/Max	£81.25
Multi 12 Channel	360	10 /chan @ Min/Mid/Max	£121.25



* A minimum of 25 pipettes are required for on-site services



Alpha Laboratories Limited
40 Parham Drive
Eastleigh, Hampshire
SO50 4NU United Kingdom

Telephone: +44 (0)23 8048 3000
Facsimile: +44 (0)23 8064 3701
Email: alpha@alphalabs.co.uk
Web: www.alphalabs.co.uk