

Certificate of Conformity

Product Code: LWXXXX

Item Description:

Lot Number: XXXXXX

Material: Low Density Polyethylene (LDPE)

Sterility: Sterile

Sterility Assurance Level (SAL): 10^{-4} (dose range 17.3 kGy to 45 kGy)

Method of Sterilisation: E-beam

Where a Sterility Assurance Level is stated, it is ensured by irradiating the product according to ANSI/AAMI/ISO 11137-2006 guidelines.

Date of Manufacture: MM - YYYY

Shelf Life: Four (4) Years

Manufacturing Location: Mexico

Other information relating to this product:

Specifications:

Length:

Bulb Draw:

Drops per ml:

Drop Volume:

Graduations:

Alpha Laboratories Ltd hereby certifies that these products have been manufactured in a controlled environment in accordance with established manufacturing guidelines. Product inspection is achieved using appropriate AQL Levels per product approved requirements. All items have been inspected and meet packaging requisites.

Certified by:



Date: 26 January 2017

Robert Vint, Managing Director